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For the construction of the vector ANep2 with the signal peptide of *A. niger* glucoamylase (ANep2_SP), ANep2_*Asc*I DNA was digested with *Nhe*I and *Asc* I and ligated with two complementary oligonucleotides, 5'-CTAGCAATGTCGTTCCGATCTCTACTCGCCCTGAGCGGCCTCGTCTGCACAGGGTTGGCAAATGTGATTTCCAGG-3' and 5'-CGCGCCTGGAAATCACATTTGCCAACCCTGTGCAGACGAGGCCGCTCAGGGCGAGTAGAGATCGGAACGACATTG-3', which encode the first 25 amino acids of glucoamylase with *Nhe*I at the 5' end and *Asc*I at the 3' end.

ASPERGILLUS MEDIA

Minimal medium(MM): for 1 liter of solid medium use:

(A) Agar	12.5g; Steam in 750 ml of demineralized or distilled water
(B) Salts	50ml of "20x Salt mix-MgSO ₄ "
MgSO ₄	0.52g of MgSO ₄ solution
C-source	10g of Glucose
TE	1ml of Trace element stock solution (see next page)
Water	

Mix (A)with (B) , dispense,supplement as needed, autoclave, and pour plates, or store at room temperature until needed.

Complete medium (CM); for 1 liter of solid CM)

To(B) of minimal medium , add the following(before agar):

Difco-Bacto Peptone	2g
Difco-Bacto yeast Extract	1g
Casamino Acids	1g
Vitamin Solution (see next page)	1ml
Water	

Mix with (A) of MM, dispense to flasks, autoclave and pour as needed.

Acetate Medium: for 1 liter of solid medium

Ammonium acetate	12g
NaCl	2g
MgSO ₄ .7H ₂ O	0.5g
KH ₂ PO ₄	3g
Trace Element	1ml
Agar	12.5g
Water	up to 1 liter

Adjust Ph to 6 by 4N-HCl ,and autoclave

Salt Mix-MgSO₄(20X stock)

To make up 1 liter, dissolve in 800 ml distilled water

NaNO ₃ (sodium nitrate)	120.0g
KCl(potassium chloride)	10.4g
KH ₂ PO ₄ (potassium monobasic)	16.3g
K ₂ HPO ₄ (phosphate dibasic)	20.9g

Hutner's Trace Element Solution : use 1 ml per liter of medium

To make up 100ml, dissolve in 80ml distilled water(in sequence, one at a time)

ZnSO ₄ .7H ₂ O(Zinc sulfate)	2.2g
H ₃ BO ₃ (boric acid)	1.1g
MnCl ₂ .4H ₂ O(Manganous chloride)	0.5g
FeSO ₄ .7H ₂ O(Ferrous sulfate)	0.5g
CoCl ₂ .6H ₂ O(cobaltous chloride)	0.16g
CuSO ₄ .5H ₂ O(Cupric sulfate)	0.16g
(NH ₄) ₆ Mo ₇ O ₂₄ .4H ₂ O(Ammonium molybdate)	0.11g
EDTA, tetrasodium salt	6.5g
EDTA, disodium salt	0.77g

Fill up to 100ml. Autoclave and store at 4°C ; unadjusted Ph will be 6.5

Vitamin Solution : use 1 ml per liter of medium

To make up 200ml dissolve in 160ml distilled water

Pyridoxin-HCl	0.02g
Thiamin-HCl	0.03g
p-aminobenzoic acid	0.15g
Nicotinic acid	0.5g
Riboflavin	0.5g
Choline-HCl	4.0g
Biotin(use conc.stock:50mg/100ml)	10ml

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